

---

# Tausand AB1000 Matlab library example: Abacus Example

## Table of Contents

Find devices and establish a connection .....	1
Read device settings .....	3
Write device settings .....	4
Read measurements from device .....	4
Multiple coincidence setting and reading .....	5
Close connection .....	6

A set of basic commands to test Tausand\_AB1000\_MatlabLibrary to be used in Matlab's command window

Author: David Guzman. Tausand Electronics, Colombia. March 2021; Last update: 15-Mar-2021.

Contact email: [dguzman@tausand.com](mailto:dguzman@tausand.com). Website: <http://www.tausand.com>

## Find devices and establish a connection

```
disp('*****');
disp('MATLAB Abacus example');
disp('*****');
disp('1. Find devices and establish a connection');

*****
MATLAB Abacus example
*****

1. Find devices and establish a connection

Find Tausand Abacus devices with reference AB1504.

disp('1a. Find Tausand Abacus devices with reference AB1504');
[ports,n] = findDevices(1504); %1504: look for only 1504 devices
disp(['Found ',int2str(n),' device(s) with reference 1504']);

1a. Find Tausand Abacus devices with reference AB1504
Progress: 1/2
Progress: 2/2
Found 1 device(s) with reference 1504

Find Tausand Abacus devices with any reference.

disp('1b. Find Tausand Abacus devices with any reference');
[ports,n] = findDevices(); %(): any device
disp(['Found ',int2str(n),' device(s) with any reference']);

if n==0
    disp('No valid devices were found. Closing.');
```

end

```
disp('Available valid devices:');  
ports
```

```
1b. Find Tausand Abacus devices with any reference  
Progress: 1/2  
Progress: 2/2  
Found 1 device(s) with any reference  
Available valid devices:
```

```
ports =
```

```
    "COM23"
```

Connect to the first device found.

```
disp('Select the first found device to connect with.');
```

```
my_device_port = ports{1}  
my_abacus = openAbacus(my_device_port)
```

```
disp('Connected to the following device:');
```

```
myidnstring = idnQuery(my_abacus);  
mydevicetype = deviceTypeQuery(my_abacus);  
disp(['idn (string): ',char(9),char(9),char(9),myidnstring]);  
disp(['device type (integer): ',char(9),int2str(mydevicetype)]);
```

Select the first found device to connect with.

```
my_device_port =
```

```
    'COM23'
```

Serial Port Object : Serial-COM23 AB1504

Communication Settings

Port:	COM23
BaudRate:	115200
Terminator:	'LF'

Communication State

Status:	open
RecordStatus:	off

Read/Write State

TransferStatus:	idle
BytesAvailable:	0
ValuesReceived:	21
ValuesSent:	7

Connected to the following device:

```
idn (string):    Tausand Abacus AB1504
```

*device type (integer): 1504*

## Read device settings

Example of reading all device settings

```
disp("*****")
disp("2. Read device settings")
disp('Settings read from device, using queryAllSettings function:');
[my_sett_data,my_sett_labels]=queryAllSettings(my_abacus);
current_settings=[my_sett_labels,my_sett_data];
disp(current_settings);
```

\*\*\*\*\*

*2. Read device settings*

*Settings read from device, using queryAllSettings function:*

<i>"sampling"</i>	<i>"2000"</i>
<i>"coincidence_window"</i>	<i>"100"</i>
<i>"delay_A"</i>	<i>"0"</i>
<i>"delay_B"</i>	<i>"20"</i>
<i>"delay_C"</i>	<i>"0"</i>
<i>"delay_D"</i>	<i>"0"</i>
<i>"sleep_A"</i>	<i>"0"</i>
<i>"sleep_B"</i>	<i>"20"</i>
<i>"sleep_C"</i>	<i>"0"</i>
<i>"sleep_D"</i>	<i>"0"</i>
<i>"config_multiple_1"</i>	<i>"224"</i>

Examples reading single settings

```
disp('Query each setting by separate:');
value=queryDelay(my_abacus,'A');
disp([' Current delay in A: ',int2str(value),'ns']);
value=queryDelay(my_abacus,'B');
disp([' Current delay in B: ',int2str(value),'ns']);
value=querySleep(my_abacus,'A');
disp([' Current sleep in A: ',int2str(value),'ns']);
value=querySleep(my_abacus,'B');
disp([' Current sleep in B: ',int2str(value),'ns']);
value=queryCoincidenceWindow(my_abacus);
disp([' Current coincidence window: ',int2str(value),'ns']);
value=querySamplingTime(my_abacus);
disp([' Current sampling: ',int2str(value),'ms']);
text=queryMultipleCoincidence(my_abacus);
if text~=""
    disp([' Current multiple coincidence setting: ',text]);
end
```

*Query each setting by separate:*

*Current delay in A: 0ns*  
*Current delay in B: 20ns*  
*Current sleep in A: 0ns*  
*Current sleep in B: 20ns*

*Current coincidence window: 100ns*  
*Current sampling: 2000ms*  
*Current multiple coincidence setting: ABC*

## Write device settings

Examples of writing a new setting value

```
disp('*****')
disp('3. Write device settings')
```

```
*****
3. Write device settings
```

Set sampling time = 2000ms = 2s

```
configureSamplingTime(my_abacus,2000);    %set sampling=2000ms
value = querySamplingTime(my_abacus);
disp([' Current sampling:',int2str(value),'ms']);
```

*Current sampling:2000ms*

Set delay in channel B = 20ns

```
configureDelay(my_abacus,'B',20);    %set delay_B=20ns
value = queryDelay(my_abacus,'B');
disp([' Current delay B:',int2str(value),'ns']);
```

*Current delay B:20ns*

Set coincidence window = 100ns

```
configureCoincidenceWindow(my_abacus,100);    %set coinc_wind=100ns
value = queryCoincidenceWindow(my_abacus);
disp([' Current coincidence window:',int2str(value),'ns']);
```

*Current coincidence window:100ns*

## Read measurements from device

```
disp('*****')
disp('4. Read measurements from device')
disp('Waiting to complete a measurement:');
waitForAcquisitionComplete(my_abacus,true);    %using default
    max_wait=10 seconds, printing messages
disp('Waiting to complete a measurement with a maxtime (1s) shorter
    than sampling time (2s):');
waitForAcquisitionComplete(my_abacus,false,1); %using
    max_wait=1second, no printing messages
```

```
*****
4. Read measurements from device
Waiting to complete a measurement:
Current ID is 183
Next data is available in    2.0s
```

```
Now, current ID is 184
Waiting to complete a measurement with a maxtime (1s) shorter than
sampling time (2s):
Warning: Maxwait expired. Consider
extending your maxwait.
```

Read data from device

```
disp('Measurements read from device, using readMeasurement
function:');
[my_meas_data,my_meas_labels]=readMeasurement(my_abacus);
current_measurements=[my_meas_labels,my_meas_data];
disp(current_measurements);
```

Measurements read from device, using readMeasurement function:

"counters_ID"	"184"
"counter_A"	"0"
"counter_B"	"0"
"counter_C"	"0"
"counter_D"	"0"
"counter_AB"	"0"
"counter_AC"	"0"
"counter_AD"	"0"
"counter_BC"	"0"
"counter_BD"	"0"
"counter_CD"	"0"
"counter_multiple_1"	"0"
"time_left"	"909"

## Multiple coincidence setting and reading

```
text=queryMultipleCoincidence(my_abacus);
if text~=""
    disp('*****')
    disp('5. Multiple coincidence setting and reading')
    %if using a device able to measure multiple coincidences
    disp([' Current multiple coincidence setting: ',text]);
    configureMultipleCoincidence(my_abacus,'ABC');
    text=queryMultipleCoincidence(my_abacus);
    disp([' New multiple coincidence setting: ',text]);
    waitForAcquisitionComplete(my_abacus);    %using default
max_wait=10 seconds
    [my_meas_data,my_meas_labels]=readMeasurement(my_abacus);
    my_index = find(my_meas_labels=='counter_multiple_1',1);
    if ~isempty(my_index)
        value = my_meas_data(my_index);
        disp([' Current coincidences in ABC: ',int2str(value)]);
    end
end

*****
5. Multiple coincidence setting and reading
Current multiple coincidence setting: ABC
```

*New multiple coincidence setting: ABC*  
*Current coincidences in ABC: 0*

## Close connection

```
closeAbacus(my_abacus)
```

*Published with MATLAB® R2017a*