
Tausand AB1000 Matlab library tester

A set of commands to test Tausand_AB1000_MatlabLibrary to be used in Matlab's command window

Author: David Guzman. Tausand Electronics, Colombia. September 2020; Last update: 15-Mar-2021

Contact email: dguzman@tausand.com. Website: <http://www.tausand.com>

```
findDevices(1504)    %searches for Tausand Abacus AB1504 devices
findDevices()        %searches for any device of the Tausand Abacus
                    family
my_ports=findDevices()    %searches for any device of the Tausand
                    Abacus family
%my_abacus=openAbacus('COM28'); %change COM port depending on your
                    results of 'findDevices'
my_abacus=openAbacus(my_ports{1}); %connects to the first found port

configureByName(my_abacus,['sampling'],[2500]) %sets sampling to 2.5s
configureChannel(my_abacus,'B',20,50) %sets delay to 20ns and sleep to
                    50ns, or the closest valid values
configureChannel(my_abacus,'a',10,59) %sets delay to 10ns and sleep to
                    59ns, or the closest valid values
configureCoincidenceWindow(my_abacus,133) %sets to 133ns, or the
                    closest valid value
configureDelay(my_abacus,'c',35) %sets to 35ms, or the closest valid
                    value
configureSleep(my_abacus,'a',35) %sets to 35ms, or the closest valid
                    value
configureMultipleCoincidence(my_abacus,"ABD")
configureSamplingTime(my_abacus,555) %sets to 555ms, or the closest
                    valid value
countersIdQuery(my_abacus)
idnQuery(my_abacus)
queryAllSettings(my_abacus)
[my_sett_data,my_sett_labels]=queryAllSettings(my_abacus)
queryCoincidenceWindow(my_abacus)
%queryDelay(my_abacus) %this line leads to an error due to a lack of
                    channel selection: a,b,c or d.
queryDelay(my_abacus,'A') %this is correctly used
queryMultipleCoincidence(my_abacus)
querySamplingTime(my_abacus)
queryAllSettings(my_abacus)
readMeasurement(my_abacus)
[my_meas_data,my_meas_labels]=readMeasurement(my_abacus)
timeLeftQuery(my_abacus)
waitForAcquisitionComplete(my_abacus)
configureSamplingTime(my_abacus,5000)
waitForAcquisitionComplete(my_abacus,true)
closeAbacus(my_abacus)
```

Progress: 1/2

Progress: 2/2

```
ans =  
  
    "COM23"  
  
Progress: 1/2  
Progress: 2/2  
  
ans =  
  
    "COM23"  
  
Progress: 1/2  
Progress: 2/2  
  
my_ports =  
  
    "COM23"  
  
  
ans =  
  
    uint32  
  
    20  
  
  
ans =  
  
    uint32  
  
    10  
  
  
ans =  
  
    uint32  
  
    130  
  
  
ans =  
  
    uint32  
  
    36  
  
  
ans =  
  
    uint32  
  
    36
```

```

ans =

    'ABD'

ans =

    uint32

    560

ans =

    uint32

    0

ans =

    'Tausand Abacus AB1504'

ans =

    11x1 uint32 column vector

    560
    130
    10
    20
    36
    0
    36
    50
    0
    0
    208

my_sett_data =

    11x1 uint32 column vector

    560
    130
    10
    20
    36
    0
    36
    50
    0

```

```

    0
    208

my_sett_labels =

    11×1 string array

    "sampling"
    "coincidence_window"
    "delay_A"
    "delay_B"
    "delay_C"
    "delay_D"
    "sleep_A"
    "sleep_B"
    "sleep_C"
    "sleep_D"
    "config_multiple_1"

ans =

    uint32

    130

ans =

    uint32

    10

ans =

    'ABD'

ans =

    uint32

    560

ans =

    11×1 uint32 column vector

    560
    130
    10

```

```
20
36
0
36
50
0
0
208
```

```
ans =
```

```
13×1 uint32 column vector
```

```
0
0
0
0
0
0
0
0
0
0
0
0
0
271
```

```
my_meas_data =
```

```
13×1 uint32 column vector
```

```
0
0
0
0
0
0
0
0
0
0
0
0
0
95
```

```
my_meas_labels =
```

```
13×1 string array
```

```
"counters_ID"
"counter_A"
```

```
"counter_B"  
"counter_C"  
"counter_D"  
"counter_AB"  
"counter_AC"  
"counter_AD"  
"counter_BC"  
"counter_BD"  
"counter_CD"  
"counter_multiple_1"  
"time_left"
```

```
ans =
```

```
uint32
```

```
49
```

```
ans =
```

```
uint32
```

```
5000
```

```
Current ID is 1  
Next data is available in 5.0s  
Now, current ID is 1
```

Published with MATLAB® R2017a